Allgemeines

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
Figure 1: Theoretical construction of a MRI <https://www.youtube.com/watch?v=0k9ILlYzmaY>

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- spin is an intrinsic form of angular momentum carried by elementary particles, composite particles (hadrons), and atomic nuclei.
- spin is like a vector quantity: it has a definite magnitude, and it has a direction.
- Spin quantum number $\frac{1}{2} \leq I$ ($I = \sum_{i=1}^A (\overline{s_i} + \overline{l_i})$). In this experiment only $I = \frac{1}{2}$, because it is easier. (^{23}Na is spin $\frac{3}{2} \rightarrow$ multiple quantum coherence)
- Number of discrete orientations (with ext. field): $2I + 1$ ($m = -I, \dots, I$)
- with $I = \frac{1}{2}$ two possible orientations ($\pm \frac{1}{2}$; spin up preferred \rightarrow lower energy) $\Rightarrow E_m = -\gamma \hbar m B_0$ γ is gyromagnetic ratio (ratio of its magnetic moment to its angular momentum)
- Protons, neutrons, and many nuclei carry nuclear spin, which gives rise to a gyromagnetic ratio $\gamma = \frac{e}{2m_p} g_n = g_n \frac{\mu_N}{\hbar}$; μ_N is the nuclear magneton, and g_n is the g-factor of the nucleon or nucleus in question. The ratio of $\frac{\gamma_n}{2\pi g_n}$, equal to $\frac{\mu_N}{\hbar} = 7622593285 \frac{\text{MHz}}{\text{T}}$
- The gyromagnetic ratio of a nucleus plays a role in nuclear magnetic resonance (NMR) and magnetic resonance imaging (MRI). These procedures rely on the fact that bulk magnetization due to nuclear spins precess in a magnetic field at a rate called the Larmor frequency, which is simply the product of the gyromagnetic ratio with the magnetic field strength. With this phenomenon, the sign of γ determines the sense (clockwise vs counterclockwise $\gamma \leq 0$) of precession. ^1H $\frac{\gamma_n}{2\pi} = 42577478518 \frac{\text{MHz}}{\text{T}}$
- Boltzmann statistic (distribution of spin up and down): $\frac{P_{\frac{1}{2}}}{P_{-\frac{1}{2}}} = \exp\left(\frac{\gamma \hbar B_0}{k_B T}\right)$; without thermal movement \rightarrow spin up; population differences results in bulk magnetisation (linearized with $k_B T$ way bigger then $\gamma \hbar m B_0$; $M = \frac{N \gamma^2 \hbar^2 I(I+1) B_0}{3 k_B T}$) parallel to the ext. mg. field

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- Zeeman effect: effect of splitting of a spectral line into several components in the presence of a static magnetic field; Zeeman sub-levels is a function of magnetic field strength
- $\Delta E = \gamma \Delta m \hbar B_0$
- NMR usually for us $\Delta m = 1$

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- with the torque the angular momentum can be reduced or raised and the rotation axis of the body can be changed
- B_0 is the external Magnetic field; the lowest energie level is, when the magnetic moments are parallel to B_0
- angular momentum: $|I| = \hbar\sqrt{I(I+1)}$
- I do not actual know what it is ment by B_0 and B_1 , but might be the coils in the experiment. B_1 coil (also gradient coil) is called the innerst coil (use: excite and detect precession) and with B_0 is probably ment the ext. field (Measurements about EFNMR $B_0 = B_E$ of the earth)
- lamor frequency $\omega_L = \gamma B_0$
- on other atoms, there is an angular momentum, which causes the precession(Lamor frequency) because of the angular momentum preservation
- Precession is the change in the orientation of the rotation axis
- if you get the item out of the angular axis(with force), the body rotates around another axis, because of the anguar momentum.
- lamor frequency $f_{Lamor} = \frac{2\gamma}{2\pi}B$; $\omega_{Lamor} = 2\pi f_{Lamor} = \gamma B = g_j \frac{q}{2m}B$
- B_1 is the magnetic field which is detected in the B_1 coil (is the most inner coil in the set up). This coil is used to detect the Signal with a LCR-circuit; B_1 is the magnetic field, which gets generated, if the magnetic mometum occure in the x-y axis

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- in conventional NMR B_0 is high so that ω_L is radio frequency (RF puls)
- also known under the word „Gradient echoes“


Figure 2: usually the spin is in the direction of the earth's field B_E or in case B_0 is applied in the direction of B_0 . In our experiment we try to manipulate this direction. We produce a RF puls so the spin flips to M_{xy} with the tip angle θ . The greyed circle shows the precession (longitudinal and transversal component; phase is random, thats why the net transversal components cancel \rightarrow net magnetisation along longitudinal axis). WE DETECT ONLY THE TRANSVERSAL COMPONENT thats why with a tip angle of $\theta = 90^\circ$ we get the best output. At $\theta = 180^\circ$ there is now output at all (under best conditions)


Figure 3: RF-Pulse with 90°


Figure 4: RF-Pulse first with 90° and then 180°

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- the bloch-equation is a macroscopic equation to calculate the nuclear magnetisation $M = (M_x, M_y, M_z)$ as a function of the relaxation time T_1 and T_2 (empirical research)(also known as equation of motion of nuclear magnetization)

- Bloch equations are for liquid and in some cases are also for solid states

$$\frac{d\vec{M}}{dt} = \gamma \vec{M} \times \vec{H}_a - \vec{e}_x \frac{M_x}{T_2} - \vec{e}_y \frac{M_y}{T_2} - \vec{e}_z \frac{M_z - M_0}{T_1} \quad (1)$$

$$\gamma \vec{M} \times \vec{H}_a = \text{is the magnetisation in the magnetic field} \quad (2)$$

- FID was first measured 1946 from Felix Bloch
- FID fourier transform nuclear magnetic resonance spectroscopy; observable NMR signal, which is generated by non-equilibrium nuclear spin; this signal gets measured in the detection coil, caused by the induction; Throughout the measuring process there is no high frequency exchange field, so this is called a „free“ degeneration; this degeneration is an exponential decay; the magnetisation can be produced by applying a short pulsed radio frequency close to the Larmor frequency of the nuclear spin; the spin is precessing, so that the induced potential in the coil is oscillating; the measured data is then often digitised and then Fourier transformed to get an NMR frequency spectrum

- the enveloping $S(t)$ from the amplitudes is described in the following

$$S(t) = S(0) \exp\left(\frac{-t}{T_2}\right) \quad (3)$$

- often the amplitude decays much faster, because of T_2^*

8.) Longitudinale Spin-Gitter- sowie transversale Spin-Spin-Relaxation, T_1 [3.3], T_2 , T_2^* [4.3]

- measure the time until the equilibrium state (GG-Lage)
- measure in 2 directions, T_1 & T_2
- T_1 : longitudinal; parallel to B_0 (in the z-Axis): time until the thermal equilibrium state (spin lattice relaxation); describes from the magnetisation $M_{x,y}$ to the M_0 state and as a result an longitudinal magnetisation M_z ; T_1 depends on the current state and the thermal equilibrium state, molecular dynamics and the strength of the magnetic field; exponentially process; the time when 63% of the magnetisation reached M_0 is called as the T_1 -relaxation time


Figure 6: area under spectrum (tells the amount of MAGNETISATION after prepolarising pulse) against the DURATION of polarising pulse. Plateau value exists. single exp. fit \rightarrow decay constant (T_1)

$$M_z(t) = M_p(1 - \exp(-t/T_1)) \quad (4)$$

$$\text{in the earth magnetic field} \quad (5)$$

$$M_z(0) \text{ equal to } M_p \quad (6)$$

$$M_z(t) = M_p \exp(-t/T_1) \quad (7)$$

- T_2 : transverse (transversale Spin-Spin-Relaxation); perpendicular to B_0 (in the x-y plane)
- T_2 is one factor of the phase coherence loss, spin-spin relaxation
- spin-spin-relaxation: caused by magnetic dipole coupling between the two neighbour-spins (in homogeneous fields, this is the main reason for phase coherence loss)

$$S(t) = S_0 \exp(-t/T_2) \quad (8)$$

$$S_0 = \text{initial magnitude at } t = 0 \quad (9)$$

for water the T_2 time constant is in the order of couple seconds


Figure 7: each echo is fourier transformed. area under spectrum (tells the size of each individual echo) at different echo times against the time. single exp. fit \rightarrow decay constant (T_2)

- in reality, the phase coherence losses come from spin-spin relaxation and the homogeneity of the magnetic field. Local homogeneity \Rightarrow a range of Larmor frequency. each nucleus has to the position associated frequency; the overall phase coherence depends of all nuclei precessing at the same frequency; frequency dispersion \Rightarrow loss of phase coherence and hence signal decay (weakens)
- $T_2^* = T_2$ relaxation + field inhomogeneities

$$\frac{1}{T_2^*} = \frac{1}{T_2} + \gamma \Delta B_0 \quad (10)$$

$$T_2^* = \text{combined relaxation effect time} \quad (11)$$

$$S(t) = S_0 \exp(-t/T_2^{ast}) \quad (12)$$


Figure 8: Visualisation of the T_1 and T_2 equilibrium state

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- spin echo trick: produce inhomogeneous field \rightarrow different ω_L . At certain time flip spins over (turn magnetisation vectors (slow spins at the front and fast spins in the back)) using transvers oscillating magnetic field \rightarrow at the startpoint they are at the same place again. ; start with 90° pulse and wait till the magnetisation decay away (spins are out of phase again) then apply 180° pulse


Figure 9: spin echo trick: point 1 (1ms) 90° pulse. Fast decay due to inhomogeneous magnetic field. At point 2 180° pulse (this is the point when the spin flips and get in sync again). at point 3 the spins (slow and fast) have different precession frequencies. At point 4 (echo) both spins are back in step. At point 5 the fast ones are at the front and the slow in the back \rightarrow out of phase again.


Figure 10: (Carr-Purcell) CPMG: This spin echo trick can be produced several times (more 180° pulses (Gaps: because have to turn off the receiver)). Amplitude dies away with time (T2 relaxation; irreversible effect); T1 is the relaxation into thermal equilibrium

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
Figure 1-5. Circuit diagram for the B_1 coil in receive mode.

Figure 11: The structural setup for the measuring method

- A LCR-circuit is the measure methode in the inner cooil of the set up
- B_1 coil an inductor with inductance L , and internal resistance R ; connected to a fixed capacitance C in spectrometer
- resonate frequency $\omega_0 = \frac{1}{\sqrt{LC}}$
- V_s is the induced voltage through the B_1 magnetic field, because of the precessing nuclear magnetism of the sample;

$$Z_{eq} = Z_R + Z_L + Z_C = R + i(\omega L - \frac{1}{\omega C}) \quad (13)$$

$$I = \frac{V_s}{Z_{eq}} = \frac{V_s}{R + i(\omega L - \frac{1}{\omega C})} \quad (14)$$

$$\xrightarrow{\omega_0 = \frac{1}{\sqrt{LC}}} I = \frac{V_s}{RS} \quad (15)$$

therefore the output voltage V_{out} at resonance of ω_0 ist

$$|V_{out}| = |IZ_C| = \frac{|V_s|}{R\omega_0 C} = \frac{|V_s|\omega_0 L}{R} = |V_s|Q \quad (16)$$

where Q is the quality factor \Rightarrow the LCR circuit amplifies the input/signal voltage

with the factor Q .

$$Q = \frac{f_0}{\Delta f} = \frac{\omega_0}{\Delta \omega} \quad (17)$$

$$= \frac{\omega_0 L}{R} \quad (18)$$

Q factor describes the width of the resonance peak

higher $Q \rightarrow$ narrower the peak \rightarrow bigger the amplification of the LCR

- in the NMR the signal to ratio (Verhältnis vom Signal zum Rauschen) is often more interesting than the signal level. The reason for this is, that beside the fact, that the signal is getting tuned by Q , also the noise gets tuned by the LCR. (Only the case by resonance).

At other frequency there is a attenuation and therefore the signal to noise ratio is boosted, of the detected signal (important: detection bandwidth is bigger than linewidth NMR)

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
Figure 1-2. A view of the three EFNMR/MRI coils

Figure 12: This is the set up for the experiment. There are three solenoids, but in the experiment they are in the biggest coil with the sample

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- 2 coaxial solenoids (Magnetspulen) and a third tube, which are 4 coils to generate a x,y,z gradients. (there are 2 x-coils, one for imaging and shimming; one for pulsed gradient NMR)
- the biggest coil has got a high current; which provides a B_p polarising magnetic field, where the nuclear magnetisation from the sample is established (in the experiment, this magnetisation is 350 times bigger than the earth magnetic field, so the signal is also 350 times bigger than in the earth magnetic field) It's activated for a few seconds, so the material gets a magnetisation and after turning off, there is still a magnetisation, which can be measured; the middle/gradient coil makes a linear magnetic field gradient (two important factors: 1.) correct the earth magnetic field by applying a small current (inhomogeneity occurs by metal in the room...), 2.) by varying the magnetic field with the coil, it is possible to produce a Larmor frequency, which depends from the γ of the nuclei; the inner coil is used to excite and detect the precession magnetisation
- shimming: adjust current in gradient coil to improve homogeneity of magnetic field (narrower peak in the spectrum)

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- earth magnetic field 50 mT \Rightarrow leads to a Larmor frequency of a few kilo Hertz for the Hydrogen atom

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
Figure 13: Additional el. mg. field to B_0 , therefore $\omega_L = \gamma B_0 + \gamma G z$. G units in $\frac{T}{m}$. $\gamma G z$ angular frequency. In 3 dimensions $G = \vec{G} = \Delta B$ and $z = \vec{r}$; simply write all with vectors


Figure 14: MRI k-space. λ gets smaller by time and k gets bigger. $\rho(z)$ says how many nuclei are at position z . $S(k)$ is the signal. 3 dimensions simply write vectors everywhere


Figure 15: Relation between $S(k)$ and the Fourier transformed (MRI image for 2 water-tubes).


Figure 16: Principle of MRI. vary G_y so that very point at the grid is scanned. After Fourier transform image is visible.

!Question: why do we need the inversion? Couldn't we apply negative G_y and G_x ??? $S(k)$ Signal looks like diffraction pattern, only after Fourier transformation the image appears.

17.) J-Kopplung; Strukturanalyse; Spektrum von 2,2,2-Trifluorethanol (CAS 75-89-8) [12.3]

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setup earths field NMR & to acquire 1H nuclear magnetic resonance (NMR) \hookrightarrow free induction decay (FID)

Terminology - MRI EFNMR: 2 coaxial solenoid + 1 solenoid (4 coils: x, y, z imaging & shimming)

B_0 (polarising) \hookrightarrow provides initial polarising field \rightarrow large current \rightarrow nuclear magnetisation

gradient coils \hookrightarrow linear magnetic field gradients

B_1 coil \hookrightarrow excite & detect precessing magnetisation ($\perp B_0$)

EFNMR controlled by digital signal processor (DSP) \cdot PC \rightarrow PROSP \rightarrow DSP

coils \rightarrow Spectrometer \rightarrow PC

multiple quantum coherence (^{23}Na (spin $\frac{3}{2}$); ^{17}O (spin $\frac{5}{2}$))

nuclei with spin quantum number $I \geq \frac{1}{2}$ possess magnetic field B_0

Number of discrete orientations (with ext. field) $\cdot 2I+1$ $m = -I, \dots, I$

This experiment $I = \frac{1}{2} \Rightarrow 2$ possible orientations (spin up $m = \frac{1}{2}$ & spin down $m = -\frac{1}{2}$)

$\Rightarrow E_m = -\gamma \hbar m B_0 \Rightarrow$ spin up preferred low energy state

\Rightarrow Boltzmann statistic: $\frac{P_{+1/2}}{P_{-1/2}} = \exp\left(\frac{\gamma \hbar B_0}{kT}\right)$ without T motions \rightarrow spin-up

\hookrightarrow population difference \rightarrow bulk magnetisation \parallel ext. mag. field

$$M = N \gamma \hbar \frac{\sum_{m=-I}^I m \exp\left(\frac{\gamma \hbar B_0 m}{kT}\right)}{\sum_{m=-I}^I \exp\left(\frac{\gamma \hbar B_0 m}{kT}\right)}$$

linear $= \frac{N \gamma^2 \hbar^2 \pm (I+1) B_0}{3kT}$ \leftarrow $\frac{H}{kT} \approx \frac{\gamma \hbar B_0}{kT} \gg 1$

Spin precession (like a spinning top) \rightarrow longitudinal & transverse component

\hookrightarrow phase is random \Rightarrow net transverse components cancel

\Rightarrow net magnetisation along longitudinal axis

NMR: 1. bulk magnetisation vector is manipulated (alternating el. mag. field, pulse)

2. transitions between energy levels ($m = \pm \frac{1}{2}$) $\Rightarrow \Delta E = \gamma \hbar B_0$

\Rightarrow Larmor frequency $\omega = \gamma B_0$

\bullet rotational NMR $B_0 \uparrow \rightarrow \omega$ is Radio frequency (RF pulse)

Earth's field NMR: ultra low $\omega \Rightarrow$ (ULF) 300-30kHz

EFNMR: ext. field is $B_E \parallel z$

B_1 RF pulse disturbs $M_z \rightarrow$ tip angle $\theta \rightarrow$ excited mag. vector along precesses with ω

B_1 detects precessing magnetisation \rightarrow free induction decay (FID)

(transverse component)

Figure 17