Measureplan for the experiment EFNMR 1 Remote

Experiment planned on July 9th 2020 at 18 o'clock

executor: Marc Neumann and Philipp Gebauer and Simon Keegan

1. login & get used to the program;
making the data file with the names (until 18:30) ☐
2. everything is connected and aligned in the earth magnetic Field ✓
3. EFNMR → **MonitorNoise** (dependent of the location and the orientation of the probe); ask tutor if the program changes C automatically; $\gamma_{H1} = 42,577 \frac{\text{MHz}}{\text{T}} \cdot 2\pi$
noise level less than $10 \frac{\mu\text{V}}{\text{Hz}}$ is okay; less than $5 \frac{\mu\text{V}}{\text{Hz}}$ is good and fewer than $3 \frac{\mu\text{V}}{\text{Hz}}$ is great (about 30 min) ☐
4. To investigate the B_1 transmit/receive coil; EFNMR → AnalyseCoil → click Analyse (note the values from the CLI);
measure the resonance frequency dependent from the Capacity;
In the figure of the script from 0 nF to 20 nF in 1/2 nF steps? (about 60 min)
5. **detect the hydrogen-signal in the water probe**
EFNMR → PulsAndCollect → measure the spektrum and change the values from B_1, C and the „shimming“ to get an better Signal; B_1 minimum and step size should be the half of the Lamor frequency (60 min) **every change and step should be noted, so you can reproduce the simulation every time**
6. **measure the longitudinal spin relaxation in the polarised magnetic field and in the earth magnetic field**
measure once τ_p in the polarized field and then t the time between the polarisation and the pulse in the magnetic field
(repeat a few times to get a good signal; 60 min) ☐
7. **measure the amplitude and the integral of the spectral peak by different shimming;** measure T_2, T_2^*
(30 min) ☐
8. **measure the puls from many spin-echos**
try to change the the time between the pulses(45 min) ☐

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